



Diameter (mm)	10	20	30	40	50	60	90	120	175	200	250	300	350	400	450
Factor	0,13	0,55	1,00	1,50	2,50	4,30	6,80	9,70	15,00	16,30	17,40	25,00	26,70	30,60	32,50

Attention!  
 Standard max. pressure load **0,1 MPa**!  
 Diagrams are only of approximate value!

Dry Air - Filterdisc 30 mm (at 20 deg. C)  
 depending on Porosity Classification



Water - Filterdisc 30 mm (at 20 deg. C)  
 depending on Porosity Classification



ROBU® Glasfilter-Geräte GmbH  
 Schützenstrasse 13  
 D - 57644 Hattert  
 Germany

Tel. : .. 49 - (0)2662 - 8004-0  
 Fax : .. 49 - (0)2662 - 8004-40  
 Email: info@robuglas.com  
 Web : <http://www.robuglas.com>