


Flowrate / Pressure Calculation

Water - Filterdisc 30 mm (20 deg. C)
depending on Porosity Classification



Conversion Factor


Diameter (mm)	10	20	30	40	50	60	90	120	175	200	250
Factor	0,13	0,55	1,00	1,50	2,50	4,30	6,80	9,70	15,00	16,30	17,40
Diameter (mm)	300	350	400	450							
Factor	25,00	26,70	30,60	32,50							

Attention !
Standard max. pressure load **0,1 Mpa !**
Diagrams only of approximate value !



Flowrate / Pressure Calculation

Dry air - Filterdisc 30 mm (20 deg. C)
depending on Porosity Classification



Conversion Factor

Diameter (mm)	10	20	30	40	50	60	90	120	175	200	250
Factor	0,13	0,55	1,00	1,50	2,50	4,30	6,80	9,70	15,00	16,30	17,40
Diameter (mm)	300	350	400	450							
Factor	25,00	26,70	30,60	32,50							

Attention !
Standard max. pressure load **0,1 Mpa** !
Diagrams only of approximate value !



ROBU Glasfilter-Geräte GmbH
Schützenstrasse 13
D - 57644 Hattert
Germany

Tel. : .. 49 - (0)2662 - 8004-0
Fax : .. 49 - (0)2662 - 8004-40
Email : info@robuglas.com
Web : http://www.robuglas.com

