



GL14 thread options




Tube, electrode,
thermocouple, etc.
6 - 7 mm diameter


GL14 thread with
open cap and
sealing ring




GL14 thread with open
cap and hose barb



GL14 thread with
open cap and
septum, for wire or
needle piercing



GL14 thread with
insert cap and
Teflon insert,
1.6, 3, or 3.2 mm ID



GL14 thread with
closed cap

Adams & Chittenden Scientific Glass
2741 Eighth Street • Berkeley CA 94710
Vox 510-843-5277 • Fax 510-843-5290
info@adamschittenden.com • www.adamschittenden.com